



LIFE  CLEAN



Mini guide

Stables, Veterinary clinics & Zoos

- **Gentler on the environment**
- **Proven effect**
- **Ready-mixed**
- **Patented solution**

For easy and efficient disinfection and cleaning of surfaces in e.g. stables, horse transports, huts, cat boxes etc. Also works for cleaning tools.

Can be used to disinfect for example:

- Stables
- Horse transports
- Huts
- Cat boxes

The product is ready-mixed and consists of a diluted aqueous solution based on ClO₂. The ready-mixed solution is not harmful to either user or nature, and can therefore be rinsed down with ordinary waste water.

To consider:

When using larger amount of LifeClean Disinfection and have poor ventilation. Remember to open windows/ventilate and, if necessary, use full or half mask with filter B (gray for inorganic gases and vapors). When soaking, use container with lid. Some sensitive people and people with asthma may experience discomfort such as coughing in poor ventilation. See more in the safety data sheet.



S. equi | EHV-1 | Equine influenza | Black mold | Parvovirus | FPLV | Coccidia spp.

For complete list see [Microbiological Efficacy Summary](#).



LIFECLEAN

AGRICULTURE

www.lifeclean.se

LIFE CLEAN

AGRICULTURE

How to apply



1 Clean away visible dirt and rinse with water (warm water is recommended).



2 Put LifeClean well on a clean cloth and wipe the surface or spray directly onto the entire surface.



3 Let work for at least 5 minutes. If necessary, rinse and dry the surface (eg. on stainless steel, leather, helmet).



Tip! When soaking tools, **note** that you should use a bucket with a lid over the solution. Let stand for 5 minutes, pick up and rinse with water. Dry.

The liquid is active as long as it is yellow. Unopened bottle in original box is valid 9 months at room temperature and 12 months in refrigerator. Store cool and avoid direct sunlight, do not freeze. Keep out of reach of children.



The bottle can be recycled

For information regarding personal protection, stability, and material compatibility please look at the actual safety data sheet.

Application techniques



Spray (Electrostatic, Pressure sprayer etc.)



LifeClean bottle & Microfiber cloth



LifeClean bottle with trigger spray

Tip!

LifeClean can be diluted up to 50% with water in regular cleaning routine.

Note!
The product should **NOT** be diluted in outbreaks such as e.g. Strangles.

